

GENADYNE
NPWT

CreedMed

UNO+

Power to Heal in the Palm of your Hand
Sustainable NPWT



- **Small and light weight** - for home or hospital use
- **Powerful** - able to address most types of wounds
- **Reusable** - multi-patient use (up to 1 year service interval)
- **Rechargeable** - No need for batteries – up to 40 hours battery usage
- **High fluid handling capacity** (200cc or 300cc canisters)
- **Choice of dressings** - supplied separately (or use regular NPWT dressings)
- **Changeable pressure settings** -40, -60, -80, -100, -125 mm HG
- **Continuous or variable** pressure settings
- **User friendly** - easy visible display with patient-safe buttons
- **Safe** - 4 audible and visual alarms (Leak, blockage, battery, and canister)
- **Multiple wound** – use the Y-connector to manage more than one wound
- **Quiet** – no disturbed nights or rest time

Reduce the waste, reduce the cost,
improve outcomes in NPWT

UNO+ Pump and Accessories

Description	Order Code
UNO+ Pump	UP-S01
200ml Canister	UP-S00D2
300ml Canister	UP-S00D3
Uno Hybrid Dressing 10x30	UP-D10130
Uno Hybrid Dressing 15x20	UP-D1520
Green Foam Dressing Kit Small	XF-DSMF1
Green Foam Dressing Kit Medium	XF-DMDF1
Green Thin Foam Dressing Kit Small	XF-DSMF1-T
Green Thin Foam Dressing Kit Medium	XF-DMDF1-T
Silver Foam Dressing Kit (antimicrobial) Small	AG-SFOAM2
Silver Foam Dressing Kit (antimicrobial) Medium	AG-MFOAM3
PVA Foam Dressing Kit Small	PVA-SFOAM1
PVA Foam Dressing Kit Medium	PVA-FOAM1
Gauze Foam Dressing Kit Small	A4-D00SP
Gauze Foam Dressing Kit Medium	A4-D00MP
Spiral Foam Dressing Kit Small	XF-SPSK1
Spiral Foam Dressing Kit Medium	XF-SPMK1



For orders please email: sales@creedmed.com

Creed Medical Ltd, Unit 2, Hamilton Court, Oakham Business Park,
Mansfield, Nottinghamshire NG18 5FB, UK

+44 (0) 1623 391578 info@creedmed.com www.creedmed.com

AWR0011/v1